



# **H<sub>2</sub>O<sub>2</sub>Bio-Decontamination      Process and      Aseptic      Transfer      of      Heat Sensitive Materials**

Sergio Mauri,  
Director, Global Marketing and Business Intelligence

**FEDEGARI**  
*GROUP*



# Presentation Outline


- Surface bio-decontamination drivers and available technologies;
- Surface bio-decontamination drivers;
- QbD, risk assessment & mitigation overview;
- Surface chemical bio-decontamination cycle;
- Cycle validation approach;
- Conclusions.





# Surface bio-decontamination drivers

- Plastic Single Use devices more and more used for Aseptic Processing → less cleaning and downtime;
- Primary containers, stoppers, clean room garments etc. outsourcing;



- Material transfer of heat sensitive items might magnify the risk of viable contamination delivery into the sterile suite;
- HPV surface bio-decontamination can mitigate this risk.




# Available technologies

- Nothing plus interlocked doors
- Spray & wipe & pray
- LAF and UV irradiation
- HPV transfer hatch
- Dry mist
- CD gas
- ND gas





# Material transfer routing



# Load patterns examples





# Risk assessment approach

What are the risks involved in an uncontrolled material transfer process?

- Health damage by supplying a contaminated product;
- Batch product rejection;
- Product shortage;
- Increased production downtime;
- Inspection failure.

These are the risks not having a good system



# Risk in aseptic processing



**Product/  
Material**

**IN/OUT**

Suite Design,  
RABS/Isolator

Equipment Design

Qualification,  
Validation

Cleaning,  
Maintenance

HVAC

Room Pressure

Adjacent Areas

Air Supply

Airflow Patterns

Personnel  
IN/OUT

SOPs

Aseptic Technique Training

Sterilization


Environment

Product

Equipment


Surface Decontamination

**FEDEGARI**  
GROUP







# Case study: stopper bag transfer





# Risk Mitigation







# What is an effective decon?

- Homogeneous distribution of the biocide within the enclosure and the load;
- Thorough and complete penetration;
- Validated biocide contact time to get the required Log Reduction of the bioburden concentration.

Any decontamination method requires a complete and thorough distribution of the biocide to get an effective decontamination.




# Decon MAL key features

- Equipment air tightness – leak test
- Turbulence HPV/Air mixture distribution
- Biocide concentration monitoring
- Equipment Chamber aeration through HEPA Filter
- Chamber positive differential pressure (15 Pa aeration and stand-by phases)
- Monitoring of Biocide effective removal below TLV requirements
- Robustness of the unit control system



# Most used biocides

- Hydrogen Peroxide vapours/mist, HPV;
- Proxy Acetic Acid dry mist, PAA;
- Chlorine Dioxide gas;
- Nitrogen Dioxide gas.





# Decontamination cycle


## A 4 phases decontamination cycle

1. Chamber pre-conditioning
2. Load of the Biocide
3. Holding time
4. Aeration






# Decontamination cycle



# Unidirectional airflow integrated into a Decontamination MAL(\*)



# MAL units examples



Surface bio-  
decontamination  
transfer hatch for  
heat sensitive  
products to be  
delivered into a Class  
A sterile suite from a  
Class C area.



**FEDEGARI**  
GROUP



# Key Factors

The right recipe to achieve a process robustness

The ingredients of our recipe are the followings:

- Material and chamber surfaces clean and dry;
- Relative Humidity;
- Temperature;
- $\text{H}_2\text{O}_2$  concentration;
- $\text{H}_2\text{O}_2$  distribution uniformity;
- Process lethality and holding time;
- $\text{H}_2\text{O}_2$  absorption/desorption;
- Aeration.




**FEDEGARI**  
GROUP

# Diagramma Psicrometrico


Unità SI

Temperatura da 0°C a +50°C

Pressione 101,325 kPa




# Parameters



TEMPERATURE/ RELATIVE HUMIDITY

# Parameters




SOURCE: PDA Europe, Pro-conference Workshop, Daikichiro Murakami, February 11, 2013





# Temperature mapping





# Biocide concentration




Most of the  $H_2O_2$  generator available on the market are controlling the biocide concentration delivered to the chamber enclosure by weight/volume rate per time;

A new vaporiser was implemented and controlled according to the instantaneous concentration of the biocide in the vapour phase.

A **PID feedback control loop** adjust the concentration within the chamber in a very narrow range.

**Better robustness, repeatability** and **decon process design** according to different loads.

# New generation $\text{H}_2\text{O}_2$ Vaporizer



**FEDEGARI**  
GROUP



# H<sub>2</sub>O<sub>2</sub> monitoring and control: two sensors

High concentration

Low concentration






# Biocide concentration

## Getting rid of stratification


- Same conditions within chamber volume;
- Uniform material decontamination;
- Avoiding cold spots.

Magnetic coupling  
propeller



Vertical inlet

Side delivery



# Aeration



## Air make up

- **Dilution** of the  $H_2O_2$  concentration with **fresh air** to strip out the biocide from material surfaces;
- **Dilution rate** is based on **150 ÷ 300** air exchanges per hour;


## RH and T

- Dry and warm air supplies the required **energy to desorb  $H_2O_2$**  from the porous loads;
- The **desorption is temperature dependent** and overcomes of the adhesion energy of the  $H_2O_2$  molecules to the surface;
- **Below 2 ppm** it worth to **cool down** the chamber and **rise the RH** for the final aeration step **below 1ppm**. Future application with large molecules of biological products will require  $H_2O_2$  residual below 50 ppb.



# Catalytic Converter

The size and expense of the HVAC system and the energy required to condition the fresh air become limiting factors and therefore air recirculation through a catalyzer bank is often required to improve the energy sustainability.







# Adsorptive materials

- The higher the  $\text{H}_2\text{O}_2$  adsorption rate of the materials, the longer the desorption and aeration time;
- Materials with higher adsorption properties can not always be avoided like Tyvek, critical issues with wet loads, paper and nylon;
- The selection of the materials influences the cycle length to the desired residual concentration;
- HEPA filters are typically an example of adsorptive material and have a real impact in the aeration time.





# Chamber tightness

Pressure decay test  
cycle start / cycle end


Dual Inflatable gaskets with  
integral integrity testing



DEFLATED



INFLATED



**FEDEGARI**  
GROUP



# Dry Cycle vs. Wet Cycle

What is the difference?

Even  $\text{H}_2\text{O}_2$  distribution  
in the enclosure;  
Lower aeration time.



DRY CYCLE?

Micro condensation – higher  
 $\text{H}_2\text{O}_2$  concentration in the liquid  
phase and shorter D values;  
RH close to the dew point.



WET CYCLE?

# Dry Cycle vs. Wet Cycle

## Pros & Cons

### Wet cycle pros:

Higher concentration of  $\text{H}_2\text{O}_2$  in the liquid phase, lower D-value, lower contact time;

### Wet cycle cons:

Longer aeration time due to the “bounce back” of  $\text{H}_2\text{O}_2$  concentration in the vapor phase;

Material compatibility to higher  $\text{H}_2\text{O}_2$  concentration in the liquid phase;

### Dry cycle pros:

More repeatable and robust process;


Even contact of  $\text{H}_2\text{O}_2$  with all exposed surfaces;

Lower aeration time;

### Dry cycle cons:

Sensible to surrounding environment T/RH variation;

Vapor phase stratification within the enclosure.



Source: Pharmaceutical Engineering  
2007

# VALIDATION PACKAGE



## Enclosure evaluation

- Load pattern study
- Air distribution study
- Enclosure pressure decay test
- Enclosure T/RH mapping
- HEPA filters leak test
- Particle counting



## Biocide uniformity evaluation

- Enclosure T mapping with biocide
- Biocide stratification
- CIs worst location



## Biological system evaluation & Optimization

- 3 BI's per location
- Enclosure D-value calculation
- **BI' worst location identification**



## Biological Performance Qualification

- 3 time BI's exposure per worst location



## Cycle safety evaluation

- **Aeration time study**
- Biocide TLV analysis





# Conclusions

- Nowadays heat sensitive disposable materials are becoming widely used in aseptic manufacturing;
- Material transfer has to be design to mitigate the risk of microbial contamination.
- Quality by Design and risk assessment approach has to be applied;
- Traditional Material Air Lock cannot guarantee a safe and validated transfer into an aseptic suite and cannot be accepted anymore;
- Surface chemical bio-contamination is now the front edge of the transfer technology;



# Thank you for your attention

And now ready for your questions...



**FEDEGARI**  
*GROUP*

Sergio Mauri  
Director, Global Marketing and  
Business Intelligence  
[SEM@fedegari.com](mailto:SEM@fedegari.com)