



ASEPTIC PROCESSING ANNEX 1 HIGHLY POTENT BIOPHARMACEUTICALS

ISPE Moscow Conference
16th and 17th of October 2017

Richard Denk
Head Sales Containment
SKAN AG richard.denk@skan.ch

Agenda

- 1 Updates on Annex 1
- 2 Highly Potent Aseptic Containment Requirements
- 3 ADC Processing. Technical Solutions for operator and product protection
- 4 Occupational Hygienic Validation/ Cleaning Validation
- 5 Question and Answer

New Annex 1

- Not yet released
- Industry is waiting for the draft to provide comments
- Latest Information from the PDA Workshop in Washington DC on the 2nd and 3rd October.
- Reasons for the update. Quality Risk Management, New Technologies like Single Used Closed Systems, Disposables, Robotics.....
- Operators should not have access to Grade A as they are the highest risk of Contamination

GMP Requirements

- No humans close to grade A (ISO 5) Area

Facilities



Conventional Aseptic Processing

Highest risk of human intervention



RABS «Restricted Access Barrier System»

Reduced risk of human intervention




Isolators

Lowest risk of human intervention

New Annex 1

- No humans close to grade A (ISO 5) Area



PDA Survey

Introduced during the Annex 1 Workshop in Washington DC on the 2nd and 3rd of October

GMP Requirements

- No humans close to grade A (ISO 5) Area

Facilities



Conventional Aseptic Processing

Highest risk of human intervention

Conventional Solution

- Operator have access to critical areas
- No Barrier
- Contamination Risk on the Curtain.
- Intensive Training and Monitoring
- Technology should be replaced to better ones.

GMP Requirements

- No humans close to grade A (ISO 5) Area



RABS “ Restricted Access Barrier System”

- Operator have access to critical areas
- Barrier but doors can be opened
- Decontamination inside of the door before closing.
- Intensive Training and Monitoring
- More and more poor designed RABS on the market.

RABS «Restricted Access Barrier System»

Reduced risk of human intervention

GMP Requirements

- No humans close to grade A (ISO 5) Area

Serious FDA Warning Letter issued to European Manufacturer of sterile Drugs, Part 2


In case of serious violations of GMP requirements, the American FDA issues a Warning Letter to the company in question. The company must react to this within 15 working days and submit a corrective action plan to the FDA.

Two aspects were criticized with regard to 21 CFR 211.113: "inadequate aseptic techniques" and "mechanical faults during media fill". The inspector supported the "inadequate aseptic techniques" with a video recording of a line set-up followed by the filling. It showed the following incorrect behaviour:

- an employee handed a pen to another employee directly above the stopper bowl
- an employee was sitting on the floor during line set-up and did not change his gown afterwards
- an employee was leaning against the cleanroom wall
- an employee left the door of an RABS open for a considerable time during the filling without working in the immediate area

GMP Requirements

- No humans close to grade A (ISO 5) Area



Isolators


- Operator have no direct access to critical areas
- Validated and accepted decontamination system with H₂O₂
- Reduced Clean Room requirements outside of the Isolator (ISO 7/8 Class C/D)
- Less Gowning of the Operator
- More and more poor designed Isolators on the market. High risk to the product

Isolators

Lowest risk of human intervention

GMP Requirements

- No humans close to grade A (ISO 5) Area



PDA Survey

Introduced during the Annex 1 Workshop in Washington DC on the 2nd and 3rd of October

Key Reason for Update

Annex 1 update

Key reasons for update

- Contamination control strategy
 - Linked to 3.6, 5.20, 5.21
 - Understand facility
 - Understand Equipment
 - Understand process
 - Update based on feedback

Chapter 5.20: New Requirements Based on the EMA Guideline on Setting Health-Based Exposure Limits in shared facilities based on the PDE „Permitted Daily Exposure“

Andrew Hopkins

Introduced during the Annex 1 Workshop in Washington DC on the 2nd and 3rd of October




Chapter 5.21: “Depending on the contamination risk, verification of cleaning of non- product contact surfaces and monitoring of air within the manufacturing area [...] in order to demonstrate effectiveness of control measures against airborne contamination or contamination by mechanical transfer.”

Agenda


- 1 Updates on Annex 1
- 2 Highly Potent Aseptic Containment Requirements
- 3 ADC Processing. Technical Solutions for operator and product protection
- 4 Occupational Hygienic Validation/ Cleaning Validation
- 5 Question and Answer

Highly Potent Biological Anti Cancer Treatments


ADCs (Antibody Drug Conj.)



Regenerative Medicine Immune-Gene-Cell Therapy



HPBs (Highly Potent Bios)



What is an Antibody Drug Conjugates ADCs and why is it so high potent ?


Warhead *Often a highly hazardous substance

Payload = Warhead + Linker

Example: Antibody Drug Conjugates



Primary mechanism of action of ADCs:
targeted delivery of a cytotoxic agent



Highly Potent Biological Anti Cancer Treatments

Protect the Patient & Protect the Operator
The overall driver for equipment & facility is the level of potency of the API payload and the level of potency of the combined product.

Operator & patient exposure needs to be understood and controlled appropriately.






Highly Potent Aseptic Products Containment Requirements

Highly Potent Aseptic Containment Requirements

The requirements placed on closed systems / containment are increasing.



Systems must be classified by means of approved limit values.



ADCs Containment Requirements

What does OEL/OEB mean?

»OEL (Occupational Exposure Limit)

Defines an average concentration load of a drug or API measured over a particular time.

The measurement is carried out in the employee's breathing area over a period of eight hours (40 hour week). The term OEL comes from the pharmaceutical industry, where internal occupational exposure limits have been calculated for a long time without being regulated by the authorities.

»OEB (Occupational Exposure Band)

It considers the toxicology of the pure substance. The aim is to provide a system categorisation that can be used to select a suitable production facility and working procedure for a product.



Containment: New Requirements Based on the EMA Guideline on Setting Health-Based Exposure Limits

New EMA Requirements

Guideline on setting health based exposure limits for use in risk identification in the manufacture of different medicinal products in shared facilities.

- Why new Requirements like the PDE (Permitted Daily Exposure) needed?
- No clear definition on which products can be manufactured in shared facilities
- Parameters like the 10ppm or 1/1000 therapeutic dose were not sufficient anymore.

New EMA Requirements

Dedicated facilities are required for manufacturing when a medicinal product presents a risk because:

- The risk cannot be adequately controlled by operational and/or technical measures,
- Scientific data from the toxicological evaluation does not support a controllable risk (e.g. allergenic potential from highly sensitising material such as beta lactams)
- Relevant residue limits, derived from the toxicologic evaluation, cannot be satisfactorily determined by a validated analytical method

Agenda

- 1 ADC Containment Requirements
- 2 ADC Processing. Technical Solutions for operator and product protection
- 3 Occupational Hygienic Validation/ Cleaning Validation
- 4 Question and Answer

ADC Process: Dispensing, Conjugation and Purification


Non Aseptic and
Toxic Bulk

Scaling and
dispensing

Conjugation

Purification

Freezing




Isolator for Purification.





Filter (Containment)
FiPa



RTP
Material in/out etc.



ADC Process: Sterile Filtration and Fill & Finish



Isolator for sterile filtration




参考画像



SART connector




ADC Process: Sterility Test



SARA
Tools/Waste etc.



ADC Process: Storage, Inspection and Packaging




ADC Process: Cleaning of Process Parts

Aseptic ADC



Off Line
Cleaning



Isolator for Cleaning

Agenda

- 1 ADC Containment Requirements
- 2 ADC Processing. Technical Solutions for operator and product protection
- 3 Occupational Hygienic Validation/ Cleaning Validation
- 4 Question and Answer






Occupational Hygienic Validation/ Cleaning Validation

How to measure Containment ?

SMEPAC Good Practise Guide


SMEPAC (Standardized Measurement of Equipment Particulate Containment)



How to measure Containment ?

Challenge: SMEPAC does not cover Aseptic Manufacturing

Aseptic Manufacturing is the Champions League on Containment.



New Method “Occupational Hygiene Validation”

1 Occupational Hygiene Validation on Fill & Finish Lines

- | | |
|------|--|
| 1.1 | Explanation of the filling line |
| 1.2 | PDE/OEL Requirements |
| 1.3 | Method of the Containment Performance |
| 1.4 | Surrogate Test Product |
| 1.5 | Risk Assessment |
| 1.6 | Used Containment Barrier |
| 1.7 | Location of the Air Samplers and Wipe Positions |
| 1.8 | Training and Good Housekeeping |
| 1.9 | Execution of the Occupational Hygiene Validation |
| 1.10 | Results / Deviation |

Cleaning Method

GMP product and operator protection and cleaning requirements of non-product-contact surfaces in aseptic Isolators.

Richard Denk SKAN AG, [\(CH\)](#)

Dr. Andreas Flückiger, F. Hoffmann-La Roche Ltd, [\(CH\)](#)

Hirokazu Kisaka, Takeda [\(JP\)](#)

Dr. Reinhold Maeck, Boehringer Ingelheim GmbH, [\(D\)](#)

Dr. Lars Restetzki, F. Hoffmann-La Roche Ltd, [\(CH\)](#)

Dr. Andreas Schreiner, Novartis, [\(CH\)](#)

Rico Schulze, GMP Inspectorate at Landesdirektion Sachsen, [\(D\)](#)

Validation of the cleanliness of non-product-contact surfaces has increased in popularity since the EMA proposed the following measures in order to demonstrate effective management of the cross-contamination risk (in Chapter 5.21 of Part 1 of its GMP guidelines): “**Depending on the contamination risk,** verification of cleaning of non- product contact surfaces and monitoring of air within the manufacturing area [...] in order to demonstrate effectiveness of control measures against airborne contamination or contamination by mechanical transfer.”

Cleaning Method

Steps of the Cleaning Method of non- product contact surfaces within aseptic Isolators:

- Risk identification based on the layout, Air Flow Simulation, routine operations within the isolator with gloves, Riboflavin Study.
- Cleaning requirements based on the ADE/PDE
 - Manual Cleaning
 - Semi automated cleaning
 - Fully automated cleaning
- Cleaning Method to demonstrate the effectiveness of the cleaning.
- Cleaning from less critical areas towards critical areas.
- Route of waste material



ISPE D/A/CH Affiliate: Containment Manual

(English Translation)



<https://www.ispe.org/publications/guidance-documents/topic>



Thank you!
Questions?



Connecting

Pharmaceutical

Knowledge

ispe.org |